

# seca 955

Digital chair scale  
with swivelling arm  
and foot rests ③



- Patients easy to move from bed to chair due to swivelling arm and foot rests.
- Fine 50-gram graduation and high capacity of up to 250 kg.
- Smooth running wheels.
- BMI function.
- Clear control element.
- Including carrying case for mains adapter.

**seca 955:**

## Consideration, comfort and safety.

**H**ow can one weigh bedridden patients or people with an impaired sense of balance as gently as possible? The answer is a seca chair scale as this can be positioned directly next to the patient, saving the patient – and the carer – unnecessary effort.

## Swivelling arm and foot rests for easy transfer to scale.

Consideration and comfort – these are particularly prominent aspects of this stable model with its circular steel tube frame. The armrests which can be folded back make a unimpeded transfer of the patient from chair or bed to the chair scale possible. The foot rests can be folded upwards and pushed under the seat. The lockable rear wheels ensure that it is safe to sit down and stand up.

## Effortless mobility.

The seca 955 has four smooth running pivoting wheels making use possible even in the narrowest of spaces.

## Clearly laid out control element.

Positioned on the back of the chair scale seca 955 and thus easy to operate is the display and control element with its large LCD display. All clearly laid out functions are easy to control with four large keys. Practical: the carrying case seca 471 provided for storing the mains adapter directly on the scale.



## Precise results for all cases.

As the zero point of the weighing function can change due to external influences – especially in the case of such a highly sensitive scale, it is automatically checked and aligned electronically once every second. And to ensure that the seca 955 determines the precise weight of restless patients, it is fitted with an effective damping system. An additional function: the BMI value and therefore the nutritional condition of the person weighed can be ascertained automatically at the push of a button.



*The ergonomically shaped handles let you push or pull the chair scale in any direction.*

**seca 955**

## Technical Data

- Capacity: 250 kg
- Graduation: 50 g < 150 kg > 100 g
- Dimensions (WxHxD):  
563 x 906 x 970 mm
- Weight: 24.5 kg
- Power supply:  
Mains adapter / Batteries
- Functions: TARE, HOLD,  
BMI, automatic switch-off
- Approval class: Ⓜ



*To make sitting down and standing up easier, the arm and footrests can be folded away.*

